



From Page No.

see page 16


a simple N-Vinyl crosslinkers, shown below may act as accelerators for matrix applications are shown below.



N-allyl maleimide



N-vinyl maleimide

a sample of maleimide  was given to M. Barkstad to test as an accelerator for matrix forming.

To Page No.

Witnessed & Understood by me,

James E. Baker

Date

Invented by

Date


Recorded by

Carl Swan

TITLE Idea

From Page No.

Run ofsted suggested converting N-Vinyl-maleimide to
N-Vinyl-sulfosuccinimide



To Page No.

Witnessed & Understood by me,

Date

Invented by


Recorded by

Date

Exhibit 2

From Page No. 536

Vinyl-Maleimide.SK2



In a 100 ml RB flask with magnetic stir bar & reflux condenser were placed 10.00875 g maleimide (lot # 90009887), 0.24960 g Na_2PdCl_4 & 35.5 g vinyl acetate (lot # 1022405). Stir & heat to refluxing. Refluxing started at 8:50 a.m. Boil. point of vinyl acetate = 72-73°C. At 1:30 p.m. - Rx turn to dark red with some solids. Continue refluxing to total 17 hours.

Refluxing stop at 1:50 p.m. - should be shut off 7 a.m. - Rx was still refluxing. Remove heating & let cool. Filter off Rx, remove excess of vinyl acetate on a Rotavap at $T=40^\circ\text{C}$ under air bleeding into the flask. We got ~15 g. residue in the flask. Add 45 ml Et_2O stir in IPA dry ice bath at $T=-20^\circ\text{C}$ for 30 min. Filter off solid, dry at RT under water aspirator to gave 5.0 g. yellow crystals /2706-21/

To Page No. 22

Witnessed & Understood by me,

Date

Invented by

Date


Dal lun

Recorded by

A. C. C. C. C.

From Page No.

Rx#1 similar, as Rx#3 in NMR tube, but using N-Vinyl Maleimide 2706-21.



F.W. = 123.11

104.06

225.15

50 mg

50.8 mg

91.41 mg

0.406 mmole

0.488 mmole

0.406 mmole

We couldn't prepare solution 50 mg N-Vinyl Maleimide in 10 mL H_2O . — Nb.

In a NMR tube was placed 50 mg N-Vinyl Maleimide & add solution of 516 mg NaHSO_3 in 10 mL H_2O . Vortex & heat at 50°C water bath for 10 min, almost all was dissolved, filter off through pipet filter to another NMR tube & submit for NMR.

Results see p. 25 back side.

Rx at RT very slow.

Rx#2 ^(0.00812M) 1 g N-Vinyl Maleimide 2706-21 + solution
100 g NaHSO_3 in 20 mL H_2O (0.0098M)
Shake at 55°C from 4 p.m. over weekend.

Rx had very small amount of solid; Rx was filtered off & water was removed with 2 x 20 mL CHCl_3 (at 60°C under water aspirator).

Got 1.71 g. yellowish residue (2706-26-1) or 93% from theory — theory yield 1.829 g.

Prepare 30 mg/0.7 mL H_2O for NMR (see p. 26 back side)

To Page No. 27

Witnessed & Understood by me,

Date

Invented by

Date

Deb Swan

Recorded by


J. C. Stehman

From Page No. Ref. 2706-21

E.E.

Vinyl-succinimide.SK2

D.B.



26.38 mg

999.99 mg

86.089 g/mole

99.088 g/mole

294.210 g/mole

125.125 g/mole

3.55 g (~3.8 ml)

1.0 g

0.025 g

1.25 g (theory)

41.28 mmole

10.0 mmole

0.0085 mmole

10.0 mmole (theory)

In a 25 mL RB flask, with magnetic stir bar were placed all ingredients. Stir & heat to refluxing. Refluxing from 3.30 p.m.

7.10 a.m. - cool Rx. Filter off through pipet filter & wash with 2x5 mL CH_2Cl_2 . Remove solvent on a Rotovap at 40°C under water aspirator with air bleeding in a flask. Got 1.3 g. yellow liquid. Add 4.5 mL Et_2O & stir in IPA -dry ice bath. Filter off solid, dry to gave 1.04 g. brownish solid / 2706-30f.

Prepare 30 mg / 0.75 mL D_2O for NMR / see p.29 back sidef.

Product looks good by NMR.

TLC was developed in $\text{CH}_3\text{OH}/\text{CHCl}_3 = 1/99$ / see p.298 / & $\text{CH}_3\text{OH}/\text{CHCl}_3 = 10/90$. We have one spot.

To Page No. 33

Witnessed & Understood by me,

Date

Invented by

Date

Dal Lwin

Recorded by

S. G. Itelman


TITLE Sulfo-N-Vinyl Maleimide Succinimide

Project No. TIPMO100
Book No. 2706

31

From Page No. _____

Ref. 2706-21



123.11
1.0 g
0.00812 M

104.06
1.02 g
0.0098 M

225.15
1.828 g (theory)
0.00812 M (-4-)

To 1.0 g N-Vinyl Maleimide (*2706-21) was added solution 1.2 g NaHSO₃ in 20 mL bi-H₂O, vortexed 5 min then placed at 55°C over on a Orbit Shaker & shaken from 2.15 p.m.

Prepare TLC, comparing Rx & starting material.

Filter off Rx-solution was slightly cloudy. Remove water with 2 x 20 mL CHCl₃, dry on a Rotovap at 60°C to give 1.67 g. light yellow crystals (2706-31).


Prepare 30 mg/0.75 mL bi-H₂O for NMR (see p. 30 back side).

Product is good.

500 mg was given to NFB for testing.

30 mg of 2706-31 was dissolved in 300 µL bi-H₂O. Added 6.0 mL of base solution - no precipitation.

1 mL Methanol + 5 mL of base soln - NaCl precipitate.



③, 30 mg of 2706-31 was dissolved in 300 µL bi-H₂O. Added 20 mL sat. K₂CO₃ - no precipitation.

Witnessed & Understood by me,

Dab Ivan

Date

Invented by

Recorded by

S. G. Gilmann

Date

To Page No. _____

TITLE N-Vinyl Maleimide

Project No. TI PMO100
Book No. 2706

39


From Page No. Ref. 2706-21

SSB

Vinyl-Maleimide.SK2

D 8

SSB



Act. wt
Maleimide
10.000/5
Na₂PdCl₄
249.69

86.089 g/mole

97.072 g/mole

294.210 g/mole

123.109 g/mole

35.5 g (~35 ml)

10 g

0.25 g

12.7 g (theory)

412.0 mmole

103.0 mmole

0.85 mmole

103.0 mmole (theory)

In a 100 ml RB flask with magnetic stir bar & reflux condenser were placed 10.0005 g Maleimide (lot # 90009887), 0.24969 g Na₂PdCl₄ & 35.5 g vinyl acetate (lot # 1022406). Stir & heat to refluxing. Refluxing started at 13.50 p.m. Boil. point of vinyl acetate = 72-73°C. Tot oil bath = 85°C.

7.15 a.m. (~17.5 hours of refluxing) - remove oil bath, let cool, filter off from solid, remove excess of vinyl acetate at 40°C with air, bleeding in a flask. We got ~14.5 g. residue in the flask. Add 45 ml Et₂O, stir in IPA-dry ice bath at T = -20°C for 30 min.

Filter off solid, dry at RT under water aspirator to gave 5.50 g. yellow crystals /2706-39/. Filtrate was stirred for 30 min more in IPA-dry ice bath at T = -20°C. Filter off, dry to gave 1.4 g. yellow crystals /39-/. Ether was removed to gave 3.0 g. yellow /2706-39/.

To Page No. 38

Witnessed & Understood by me,

Date

Invented by

Date

Dab Swan

Recorded by


S. S. S. S.

back side.

From p. 39.
solids (2706-398). Seems that product started to polymerize. oil (222) is 95% pure by char

Redissolve solid (39-3) in 25 mL CHCl_3 by shaking on an Orbit Shaker for 20 min, filter off solids that didn't dissolve.

Remove CHCl_3 on a Rotovap at RT under water aspirator, with air bleeding into a flask. Tracey of solvent were removed by sweeping on with air, to gave 1.41 g. yellow solid
/2706-39-38/.



666


Exhibit 10

From Page No.

see pages 16-19

Purpose: to determine if *X* formamide would be a solvent for the reaction of *N*-Vinylformamide and the potassium salt of sulfopropylacrylate.

NVF-rx2.sk2



Procedure: The ingredients were stirred at an unknown temperature (25 to 90 C must likely). After 20 hours 0.1 ml was treated with 0.5 ml methanol and 0.5 ml chloroform. Removal of the volatiles gave 99 mg residue 2683-30-1 (mainly formamide any product?). The residue was washed with a second portion of methanol 0.5 ml and chloroform 0.5 ml. The clear liquid was again removed and evaporated to give 2683-30-2 (12.9 mg). The residue after two washings was dried to give 2683-30-3 (6.4 mg). Three samples were made for NMR comparisons- potassium sulfopropylacrylate 2683-30-4, formamide 2683-30-5, and *N*-vinyl formamide 2683-30-6. A final reaction sample 0.1 ml worked up with methanol and chloroform was labeled 2683-30-7. Sample 1, 2 and 7 appeared to show a new four lined NMR peak at -6.95 ppm. This new NMR peak may be evidence for the presence of the desired product.

Job

To Page No.

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Seonul Bae

Dab Swan

SurModics Intellectual Property and Proprietary Product Idea Form

②

Originator(s)

Date

Ron Ofstead and Dale Swan

Title/Key Words

N-vinylamides as accelerators in matrix formation

Reference (Personal Notes/Notebook Number and Pages)

2683-16,20,26

Brief Description

Cells can be covered with a protective hydrogel coating. The polymerization of PEG-triacrylate around the cells is accelerated by the addition of N-vinylamides. In addition the presence of sulfonate containing monomers (ie AMPS) have been useful in improving biocompatibility. The idea was to synthesize reagents containing N-vinylamides and sulfonate functionality. The attachment of figures 1 to 4 show the reactions used to make N-vinyl amides.

Advantages and Features

The materials proposed can be made in one or two steps from available materials. Preliminary tests indicated firm gels resulted from the cyclic products synthesized.

Reduced to Practice (Date/Notebook Number and Pages)

2706-21, 26, 30, 31, 37, 39 from

Submitted by

Dale Swan

DALE SWAN

Signature

Printed Name

Originator(s)

Date

R. Ofstead

R. Ofstead

Signature

Printed Name

Originator(s)

Date

Read and Understood by

Anthony Dalluied

Anthony C. Dalluied

Signature

Printed Name

Witness

Date

Jerome C. Behrens

Jerome C. BEHRENS

Signature

Printed Name

Witness

Date

PROPRIETARY
SurModics, Inc.

Exhibit 12

Based on the two previous batches made (#1 looked great; #2 didn't work in their system)

~~Examine~~ ^{Examine} ~~and~~ ^{and} ~~on~~ ^{on} ~~server~~ ^{server} ~~experiments~~ ^{experiments} ~~to~~ ^{to} ~~test~~ ^{test} ~~the~~ ^{the} ~~system~~ ^{system} ~~of~~ ^{of} ~~the~~ ^{the} ~~accelerant~~ ^{accelerant}.
 Make up ~~the~~ ^{the} ~~solutions~~ ^{solutions} ~~@~~ [@] ~~different~~ ^{different} ~~levels~~ ^{levels}
 of sub-irradiation-succinate. Lot # 2703-K-(1,2,3,4,5)

received ~ 50mg of each. Added ~ 5mg of each to 3% HA, 0.25% MSA solution, & let
 mix for 1 hour on 37°C shaker (Ammonia vials labeled 1-5, for representational - see previous
 page for vial set-up)

Net: mixing, using 75ul to make better mix, & illuminate for 45 sec.

- 1) Soft, no matrix, bleaching
- 2) Soft, matrix, no bleaching
- 3) Soft, firm matrix,
- 4) ~
- 17) ~

5 solutions set o/N @ Room Temperature. All solutions, when illuminated, looked similar as first ones

during mixed o/N @ 37°C shaker - when illuminated, solutions looked similar; #3 may have
 been a bit softer, but hard to tell.

